

# BlisLid MA

PET lidding for blister applications




Lidding films to round off our one-stop-shop-portfolio: Liveo Research offers now also sealing films from 20 to 50 microns that allow our customers to have a One-Material Blister (OMB) when choosing one of our perfect matching basefilms.

## AVAILABLE STRUCTURES


### BlisLid MA (PL)

Extruded PET lidding film as peel-solution



### BlisLid MA (PT++)

Extruded PET lidding film as push through solution



## FIELDS OF APPLICATION

- Ethicals
- Nutraceuticals
- Supplements
- Generics
- Probiotics
- Veterinary
- OTC
- Dietary
- more upon request

## QUALITY & REGULATORY

- Current valid EU and US Pharmacopeia
- CFR 21
- EC 1935/2004, (EU) No. 10/2011
- REACH – Regulation (EC) No. 1907/2006
- US CONEG Regulations
- ISO 15378

## PROPERTIES & BENEFITS



### HALOGEN FREE

film solution without halogens included



### USER FRIENDLINESS

easy push through properties available comparable to aluminium lids



### EFFICIENCY

machine speed comparable to PVC/Alu-blister production



### PLUG & PLAY

no machine modifications required, minor changes to process conditions



Innovation meets experience – development meets customers demands: Benefit from our latest developed PET films with further improved properties! Contact us for further information!

[info@liveoresearch.com](mailto:info@liveoresearch.com)  
Follow us on LinkedIn:

[www.liveoresearch.com](http://www.liveoresearch.com)  
[www.linkedin.com/company/liveo-research/](https://www.linkedin.com/company/liveo-research/)